

Cell-Lysate Experiment of CCH2:

Materials Needed:

250 mL Erlenmeyer flasks
M9 minimal salt media
Chloramphenicol (25mg/ml)
1M IPTG
200 mM 6-CNA dissolved in methanol
250 mL Erlenmeyer flasks with 50 mL LB
Spectrophotometer
Cuvette
20mM (+80mM NaCl) MOPS Buffer
7M Urea Buffer
Sonicator
1.5mL Eppendorf tubes
2 mL Eppendorf Tubes
10mL glass test tubes
Ethyl Acetate
12M Hydrochloric Acid
pH paper
GC vials
GC vial micro-inserts
GC vial caps
Nitrogen Gas
Pyridine
BSTFA + TMSC 99:1 (Derivatizing Agent)
Vortex
Centrifuge
GC/MS
Pipet

M9 minimal salt media components:

33.9g/L Na_2HPO_4
15g/L KH_2PO_4
5g/L NH_4Cl
2.5g/L NaCl

Protocol:

- 1) See Determination of CCH2 rate protocol and repeat steps 1-6 [\[1\]](#).
- 2) Four aliquots are taken of 1.5mL and transferred into four Eppendorf tubes.
- 3) The samples were centrifuged at 15000rpm for 30 minutes at 4°C and the supernatant is discarded.
- 4) The cell pellets are resuspended in 800 μL of MOPS.
- 5) Two of the four resuspended pellets were sonicated at 1W two times. The other two were sonicated on the lowest setting four times. The aliquots

- were sonicated for ten seconds and had one minute of wait time before the subsequent sonication, while being kept on ice.
- 6) One of each sonication time samples is centrifuged at 15000rpm for 30 minutes at 4°C and the supernatant transferred to Eppendorf tubes.
 - 7) The other sonication time samples stay on ice without being centrifuged.
 - 8) 380 µL of each sample are transferred into Eppendorf tubes and 20 µL of 6-CNA was added.
 - 9) The samples are incubated overnight at 30°C.
 - 10) The supernatant fractions are extracted via the liquid-liquid extraction protocol, with modified quantities of 12M HCL and Ethyl Acetate [2].
 - 11) The samples were then prepared via the GC-MS sample preparation protocol [3].

Other Protocols:

[1] Step 1-6 of protocol- Determination of CCH2 rate:

[2] Liquid-Liquid Extraction for isolation of 6-Chloronicitinoic Acid and 6-Hydroxynicitinoic Acid

[3] Gas Chromatography-Mass spectroscopy conditions to run 6-Chloronicitinoic Acid and 6-hydroxynicitinoic Acid and sample preparation: