

Curdlan Quantification :

Overview :

Curdlan and other (1→3)- β -D-glucans specifically bind the triphenylmethane dye, aniline blue, and a benzophenone fluorochrome found in the dye mixture. Calcofluor and Congo Red also bind to curdlan and induce fluorescence, but these dyes are not specific for (1→3)- β -glucans. So we used aniline blue dye to quantify our curdlan production from *E.coli* and *S.cerevisiae*.

Procedure :

References : 1,3--Glucan Quantification by a Fluorescence Microassay and Analysis of Its Distribution in Foods ; YUAN-TIH KO* AND YU-LING LIN, J. Agric. Food Chem. 2004, 52, 3313–3318

1. Prepare DYE MIX :
 - 40 volumes of 0.1% aniline blue in water
 - 21 volumes of 1 N HCl
 - 59 volumes of 1 M glycine/NaOH buffer pH 9.5
2. Diluted sample (also Curdlan Standard) 10-50-fold with 1 N NaOH to a final volume of 300 μ L in a 1.5 mL microcentrifuge tube
3. Add 30 μ L of 6 N NaOH
4. Incubated at 80 °C for 30 min
5. The tube was immediately put on an ice bath
6. Add 630 μ L DYE MIX into the tube and mix
7. Incubated at 50 °C for 30 min
8. The unbound fluorescent dye was decolorized at room temperature for 30 min
9. Measure the fluorescence intensity with a fluorescence spectrophotometer
 - emission wavelength of 502 nm (20 nm slit width)
 - excitation wavelength of 398 nm (20 nm slit width).