

***E.coli* Curdlan Purification**

Overview :

Curdlan molecules may have as many as 12,000 glucose units and are insoluble in water, alcohols and most organic solvents, while they are soluble in a diluted alkali solution (0.25 M NaOH).

We have two protocols for Curdlan purification from *E.coli* based on this property: curdlan is soluble in a NaOH solution.

Procedure n°1 :


Reference : *Demonstration of Curdlan-type polysaccharide and some other in microorganisms with aniline blue, Itaru Nakanishi, Tokuya Harada and all, J. Gen. App!. Microbiol., 22,1-11(1976)*

- 1) Mix the culture with an equal volume of 1 N NaOH and then centrifuge at 10,000 rpm for 10min to remove the cells.
- 2) Neutralize the resulting supernatant by adding 3 N HCl
- 3) A precipitate is formed. Collect it by centrifugation at 10,000rpm for 10min and wash three times with water by centrifugation. Then dehydrate with acetone and dry it in vacuo.
- 4) This precipitated polymer is the curdlan-type polysaccharide.

Procedure n°2 :

Reference : « Exopolymers from curdlan production: incorporation of glucose-related sugars by *Agrobacterium* sp. strain ATCC 31749. Jin W. Lee & all, Canadian Journal of Microbiology, 1997, Vol. 43, N° 2 : pages 149-156 »

Fig. 1. Flow diagram depicting the purification procedure for curdlan.
SN, supernatant.
Adapted from the reference.



- 1) Centrifuge the culture at 12000 x g for 30 min.
- 2) Add the pellet to an equivalent volume of 0.5 N sodium hydroxide at 3°C
- 3) Stir the mixture for 10 min and let it stand for 3 h at the same temperature.
- 4) Centrifuge the resulting viscous solution at 12000 x g for 40 min
- 5) Precipitate Curdlan in the clear supernatant by neutralization with 10% acetic acid
- 6) Repeatedly wash with DWater, acetone and ether.