

Restriction analysis of the three plasmids, integration of dCas9_VP64

In charge: Yeast group (Cyril, Loïc, Vincent, Axel, Joseph)

Start date: 30/05 - End date: 30/05 OR Date: 01/06

Aim:

Digest the three plasmids pTPGI_dCas9_VP64, pRPR1_gRNA_handle_RPR1t and plp with restriction enzymes. Then, the electrophoresis will tell if we have the right plasmids.

Schedule :

- Digest the three plasmids with restriction enzymes. Run agarose gel electrophoresis.

Results, discussion and comments:

Step 1: Restriction analysis to check the presence of the three plasmids

Date: 30/05

Duration: 3 hours

Participants:

Axel, Cyril, Loïc

Aim:

Until now, we can't be sure that pTPGI_dCas9_VP64, pRPR1_gRNA_handle_RPR1t and plp are really present in our samples. This step allows to verify by checking the size of different fragments of the plasmid.

Material:

- NEB Buffers (2.1 and 3.1)
- REs (HindIII, XhoI, PstI, BsrBI, EcoRV)
- Plasmid DNA (dCas9 VP-64, gRNA, plp)

Protocol:

Based on [NEB protocol](#) (and LISV p62).

Master mix for each plasmid :

- Add 1/10th (1 µL) of the final volume of 10X concentration buffer (NEBuffer 3.1).
- Add water to reach a final total volume of 10 µL.
- Add 1 µg of plasmid DNA.
- Add 1 µg (1 µL) of each restriction enzyme. Restriction enzymes should be last component added to reaction. Keep the restriction enzymes on ice if they leave the freezer.

Samples were loaded according to [calc sheet T1 E1 Restriction Analysis](#)

Most enzymes were 100% efficient with NEBuffer 2.1 and 3.1. PstI was 75% efficient with NEBuffer 2.1.

- Mix carefully by tipping the tube or pipetting (no vortexing) followed by a quick spin.
- Incubate your reactions at 37°C for 60 min in heating block.
- Put samples on ice until agarose gel is ready.

Comments:

pTPGI_dCas9_VP64 : HindIII, PstI, XhoI

plp : PstI, HindIII, BsrBI

gRNA : EcoRV, XhoI and BsrBI

If restriction enzymes come from NEB, we have to use the NEB protocol with 10X NEBuffer instead of FastDigest Green Buffer.

Results:

Expected results: See PDF files "Expected Gel.pdf".

The first lane is for undigested plasmid.

Miniprep of pTGI_dCas9_VP64

Date: 21/07

Participants:

Joseph

Aim:

Isolate the pTGI_dCas9_VP64 plasmids.

Material:

- Overnight cell culture from single colony inoculation
- QIAprep Spin Miniprep Kit Components:
 - P1 = resuspension buffer + RNase A and LyseBlue
 - P2 = lysis buffer
 - N3 = neutralizing buffer
 - PB = wash buffer
 - Pe = wash buffer
 - EB = elution buffer
 - QIAprep spin column
 - Collection tube
 - 1.5 microfuge tube

Protocol:

protocol from Qiagen Miniprep Kit.

Comments:

We used 3 cultures of 3ml of overnight culture, we pooled them at the beginning of the protocole. 2 of the cultures came from one glycerol stock and one from the other.

Results:

The following concentration were obtained :

- 117.1 ng/ul. 260/280 : 1.9 260/230 : 1.68 -> There may be residual guanidine

Storage: In the -20°C freezer.

Prepare for pTPGI_dCas9_VP64 integration

In charge: Yeast group (Cyril, Loïc, Vincent, Axel, Joseph)

Start date: 15/07 - End date: --/--

Aim:

Integrate pTPGI_dCas9_VP64 into yeast genome so that the dCas9_VP64 protein is always expressed by the cells.

Schedule :

- Prepare medium and plates
- Linearize the plasmid
- Integrate in yeast

Step 1: medium preparation

Date: 22/07

Participants: Loic, Joseph

Aim:

Prepare the needed plates and medium.

Material:

Step 2: Test potential single cutters for plasmid linearization

Date: 22/07

Participants: Loic, Joseph

Aim:


restriction analysis in order to choose the best single cutter.

Material:


Procedure:

Samples loaded on two gels. Twice too many enzymes but we don't give a shit.

lane	1	2	3	4	5	6	7
enzyme (20 U)	NotI HF	NotI HF	EagI HF	EagI HF	EagI HF	KpnI HF	no
buffer (1 uL)	2.1	cutsmart	2.1	3.1	cutsmart	cutsmart	no
plasmid	pTPGI_dCas9_VP64 (100 ng)						
water	7 uL						8 uL



lane	1	2	3	4	5	6
enzyme (20 U)	no	PciI	Acc651	NaeI	StuI	StuI
buffer	no	3.1	3.1	cutsmart	2.1	cutsmart
plasmid	pTPGI_dCas9_VP64 (100 ng)					
water	8 uL	7 uL				



Results:

Step 3: linearize the plasmid

Date: 22/07

Participants: Loic, Joseph

Aim:

Check linearization.

Material:

Procedure:

Samples were loaded on the gel for verification.

lane	1	2	3
enzyme (40U=2uL)	no	EagI HF	NotI HF
buffer (5 uL)	no	cutsmart	cutsmart
plasmid pTPGI_dCas9 _VP64	4 ug	4 ug	4 ug
water	up to 50 uL		

Step 4: transformation protocol

Date: 22/07

Participants: Loic, Joseph

Aim:

integration into genome

Material:

Procedure:

EPFL igem 2014 → [protocol of Julie](#)

Result :

transformed yeast on agar plates.

Start of inoculation : 22/07 at 18h.