

Lab Notebook:  
Cellulose Project

7/6/2015- 7/10/2015

1. Grew culture of *G. hansenii* in medium
2. Still growing for the entire week

7/13/2015 – 7/17/2015

1. Culture still growing
2. Read more about cellulose binding domain

7/20/2015 – 7/24/2015

1. Designed the cellulose binding system – originally plan to do antibody design. Reason is that the antibody can bind to one specific region while chemically activated streptavidin can bind to multiple regions of the proteins and can inhibit the function of the protein.
  - a. Cons: The sequence for the Flag tag antibody is not sequenced thoroughly meaning that there are high risk the design will not work just because the sequence is wrong.

Antibody Design:

T7 - rbs - pel B - ATG - sfGFP - linker - CBD- linker - VH - linker - VL - (His)6tag - Terminator

T7 promoter - rbs - sfGFP - linker - CBD - linker - Pel B - ATG - Vh - Linker - VI - (His)6Tag - Terminator

711 bp + 700 +

ATGAAATACCTGCTGCCGACCGCTGCTGCTGGTCTGCTGCTCCTCGCTGCCCAGCCGGCGATGGCC + ATG +  
aatgctacgccaactaagggtgcaacccgaccaacacagcaacgcctacaaaaagcgctacagcaacacctacaagaccgtcagttcctacaaac  
acac  
cgactaacacaccggcaatacacctgtttcaggcaacttgaaggct gaGttc  
tataactcaaataccgagtatacaactaacagtattaatccgcagtt  
taaagtaacaaatacaggatcaagtgaattgatctttcaaagcttacattgagatactattacaccgttgatggccagaaggaccagactttctggtg  
t  
gaccatgcagctatcataggtagcaacggctcatacaacggcatcacatcaaatgtaaaaggaaacattcgtaaagatgagctcaagcacaataac  
gcag  
acacatacctcgaaataagtttcacaggtggcactttggaacctgggtgctcatgtacagatacagggtaggtttgcgaaaaatgactggagtaattata  
c  
acagtcaaatactactcatttaagtcagcatcacagttcgtagaatgggatcaggttacagcatatttgatggagtacttgatggggtaaagaacc  
a  
ggaggatcagtagttccgtcaacacagccggttaacaacccaccggcaacaaccaagccgccagcaacaacaaaccaccggctaccacgattcct  
ccat  
cagacgatccg + VH + linker + VL + (His)6Tag

taatacgaactcactataggagatactagagaaaggagagaataactagatggccggc +  
ATGAAATACCTGCTGCCGACCGCTGCTGCTGGTCTGCTGCTCCTCGCTGCCCAGCCGGCGATGGCC  
+ ATG +  
cgtaaaggcgaagagctgttactggtgtcgtccctattctggtggaactggatggtgatgtcaacggctcataagttttccgtgcgtggcgagggtgaa  
g  
gtgacgcaactaatggtaaactgacgctgaagttcatctgtactactggtaaactgccggtaccttggccgactctggtaacgacgctgacttatggtgt  
tcagtgtttgtcgttatccggaccatatgaagcagcatgacttctcaagtcgccatgccggaaggctatgtgcaggaacgcacgatttcctttaag  
gatgacggcacgtacaaaacgctgcggaagtgaatttgaaggcgataccctggtaaaccgcattgagctgaaaggcattgactttaagaagac  
ggca

atatcctgggccataagctggaatacaattttaacagccacaatgtttacatcaccgccgataaacaataatggcattaaagcgaatttttaaatt  
cg  
ccacaacgtggaggatggcagcgtgcagctggctgatcactaccagcaaaacactccaatcggtgatggctctgttctgctgccagacaatcactatc  
tg  
agcacgcaaagcgttctgtctaaagatccgaacgagaaacgcgatcatatggttctgctggagttcgtaaccgcagcgggcatcacgcatggtatgga  
tg  
aactgtacaaa +  
aatgctacgccaactaaggggtgcaaccccgaccaacacagcaacgcctacaaaaagcgctacagcaacacctacaagaccgtcagttcctacaaac  
acac  
cgactaacacaccggcaaatacactgtttcaggcaacttgaaggtcgaattctataactcaaatccgagtatacaactaacagtattaatccgcag  
tt  
taaagtaacaaatacaggatcaagtgaattgatctttcaaagcttacattgagatactattacaccgttgatggccagaaggaccagactttctggtg  
t  
gaccatgcagctatcataggttagcaacggctcatacaacggcatcacatcaaatgtaaaaggaaacattcgtaaagatgagctcaagcacaataac  
gcag  
acacatacctcgaaataagtttcacaggtggcactttggaacctgggtgctcatgtacagatacagggtaggtttgcgaaaaatgactggagtaattata  
c  
acagtcaaatgattactcatttaagtcagcatcacagttcgtagaatgggatcaggttacagcatatttgaatggagtacttgatggggtaaagaacc  
a  
ggaggatcagtagttccgtcaacacagccggttaacaacccaccggcaacaaccaagccgccagcaacaaccaaccaccggctaccacgattcct  
ccat  
cagacgatccg + VH + linker + VL + (His)6tag - Terminator

T7 promoter - rbs - sfGFP - linker - CIP - linker -

taatacgactcactatagggagatactagagaaaggagagaaatactagatggccggc  
cgtaaaggcgaagagctgttactggtgtcgtccctattctg  
gtggaactggatggtgatgtcaacggtcataagtttccgtgctggcgagggtgaaggtgacgcaactaatggtaaactgacgctgaagttcatctgt  
a  
ctactggtaaactgccggtaccttggccgactctggtaacgacgctgacttatggtgttcagtgcttctgctgttatccggaccatatgaagcagcatga  
cttcttcaagtcgccatgccggaaggctatgtgcaggaacgcacgatttctttaaggatgacggcacgtacaaaacgcgtgcggaagtgaatttg  
aa  
ggcgataccctggtaaaccgcattgagctgaaaggcattgactttaagaagacggcaatatcctgggccataagctggaatacaattttaacagcc  
aca  
atgtttacatcaccgccgataaacaataatggcattaaagcgaattttaaaattcgccacaacgtggaggatggcagcgtgcagctggctgatca  
cta  
ccagaaaacactccaatcggtgatggctctgttctgctgccagacaatcactatctgagcacgcaaagcgttctgtctaaagatccgaacgagaaac  
gc  
gatcatatggttctgctggagttcgtaaccgcagcgggcatcacgcatggtatggatgaactgtacaaa accggc  
aatgctacgccaactaaggggtgcaa  
ccccgaccaacacagcaacgcctacaaaaagcgctacagcaacacctacaagaccgtcagttcctacaaacacacccgactaacacacggcaaata  
cacc  
tgtttcaggcaacttgaaggtcgaattctataactcaaatccgagtatacaactaacagtattaatccgcagtttaagtaacaaatacaggatcaag  
t  
gcaattgatctttcaaagcttacattgagatactattacaccgttgatggccagaaggaccagactttctggtgtgaccatgcagctatcataggttagca  
acggctcataacggcatcacatcaaatgtaaaaggaaacattcgtaaagatgagctcaagcacaataacgcagacacatacctcgaaataagttt  
cac  
agggtggcactttggaacctgggtgctcatgtacagatacagggtaggtttgcgaaaaatgactggagtaattatacacagtcaaatgattactcatttaa  
g

tcagcatcacagttcgtagaatgggatcaggttacagcatatttgaatggagtacttgatggggtaaagaaccaggaggatcagtagttccgtcaac  
ac  
agccggtaacaacccaccggcaacaaccaagccgccagcaacaacaaaccaccggctacca cgattcctccatcagacgatccg

sfGFP: (711 bp)

cgtaaaggcgaagagctgttactgggtgctgctccatttctgggtggaactggatgggtgatgtcaacggtcataagttttccgtgctggcgagggtgaa  
g  
gtgacgcaactaatggtaaactgacgctgaagttcatctgtactactggtaaactgccggtaccttggccgactctggtaacgacgctgacttatgggtg  
tcagtgtttgctcggtatccggaccatatgaagcagcatgacttctcaagtcgccatgccggaaggctatgtgcaggaaacgcacgatttcctttaag  
gatgacggcagctacaaaacgctgcggaagtgaatttgaaggcgataccctggtaaaccgcattgagctgaaaggcattgactttaagaagac  
ggca  
atatctggggcataagctggaatacaattttaacagccacaatgtttacatcaccgccgataaacaacaaaaaatggcattaaagcgaattttaaaatt  
cg  
ccacaacgtggaggatggcagcgtgcagctggctgatcactaccagcaaaacactccaatcggtgatggctctgttctgctgccagacaatcactatc  
tg  
agcacgcaaagcgttctgtctaaagatccgaacgagaaacgcgatcatatggttctgctggagttcgtaacgcgacgggcatcacgcatggatgga  
tg  
aactgtacaaa

linker - Cipa - linker

aatgctacgccaactaagggtgcaaccccgaccaacacagcaacgcctacaaaaagcgctacagcaacacctacaagaccgtcagttcctacaaac  
acac  
cgactaacacaccggcaaatacacctgtttcaggcaacttgaaggtcgaattctataactcaaatccgagtgatacaactaacagtattaatccgcag  
tt  
taaagtaacaaatacaggatcaagtgaattgatctttcaagcttacattgagatactattacaccgttgatggccagaaggaccagactttctgggtg  
t  
gaccatgcagctatcataggtagcaacggctcatacaacggcatcacatcaaatgtaaaaggaaacattcgtaaagatgagctcaagcacaataac  
gcag  
acacatacctcgaaataagtttcacaggtggcactttggaacctggctcatgtacagatacagggtagggttgcgaaaaatgactggagtaattata  
c  
acagtcaaatgattactcatttaagtcagcatcacagttcgtagaatgggatcaggttacagcatatttgaatggagtacttgatggggtaaagaacc  
a  
ggaggatcagtagttccgtcaacacagccggtaacaacccaccggcaacaaccaagccgccagcaacaacaaaccaccggctaccacgattcct  
ccat  
cagacgatccg

PelB Sequence:

ATGAAATACCTGCTGCCGACCGCTGCTGCTGGTCTGCTGCTCCTCGCTGCCAGCCGGCGATGGCC

heavy chain: domain 3 1-118

evqlqqsgge lakpgasvkm sckssgyft ayaihwakqa agaglewigy iapaagaaay  
61 naafkgkatl aadkssstay maaaaltsed savvycaraa aagadywqgq tltlvssak

Linker is GGGGSGGGSGGGGS

light chain: 2-111 out of 216 aa

vlmtqaplt lpvslgdqas iscrssqaiv hangntylew ylkpgqspa lliykvanrf  
61 sgvpdrfsgs gsgtdftlki srveaedlgv yycfqgahap ytfgggkclei

domain 1: 1-114

dvlmtqapltlpvslgdqasiscrssqaivhangntylewylqkpgqspa lliykvanrf  
61 sgvpdrfsgs gsgtdftlki srveadlgv yycfqgahap ytfgggtkle ikra

His6 tag:

flag tag: DYKDDDK

anti-Flag M2

- 2. Design of the streptavidin-biotin construct:** The streptavidin construct has one of the strongest binding affinity with biotin. Further research shows that bira can allow biotinylation of a specific peptide acceptor region known commercially as the Avitag sequence. The design of the construct is detailed below:

Streptavidin-Biotin Construct:

- (1) CIP - Streptavidin - terminator
- (2) Pro - RBS - AP Tag - Protein of Interest - end
- (3) Pro - RBS - BirA - end

BIRA - in cytoplasm expressed from <http://www.uniprot.org/uniprot/P06709.fasta>  
MKDNTVPLKLIALLANGFHSGEQLGETLGMSRAAINKHQTLRDWGVDFVTPGKGYSL  
PEPIQLLNQKQILGQLDGGSVAVLPVIDSTNQYLLDRIGELKSGDACIAEYQQAGRGRRG  
RKWFSPFGANLYLSMFWRLEQGPAAAGLSLVIGIVMAEVLRLGADKVRVKWPNDLYLQ  
DRKLAGILVELTGKTGDAAQIVIGAGINMAMRRVEESVVNQGWITLQEAGINLDRNTLAA  
MLIRELRAALELFEQGLAPYLSRWEKLDNFINRPVKLIIGDKEIFGISRGIDKQGALLL  
EQDGIKPWMGGEISLRSAEK

Avitag-Can be at the N or C terminus

GLNDIFEAQKIEWHE

GLNDIFEAQKIEWHE

if fused at the N terminus, need a methionine start codon with a second codon  
consistent with efficient translation and protein stability.

starting MA sequence, ATGGCT,

85-11 tag: <http://www.ncbi.nlm.nih.gov/pmc/articles/PMC2144313/pdf/10211839.pdf>  
MAGGLNDIFEAQKIEWHE

At C-terminus, the 15-mer is connected to the protein through a minimum a GG peptide  
linker

Streptavidin : Part BBa\_J36841

gctgaagctggtatcaccggcacctggtacaaccagctgggatccacctcatcgttaccgctggtgctgacggtgctctgaccggtacctacgaatcc  
g  
ctgttggttaacgctgaaagccgctacgttctgaccggtcgttacgactccgctccggctaccgacggtccggaaccgctctgggttgaccgttgcttg  
gaaaaaactaccgtaacgctcactccgctaccacctggtctggccagtagcttggtggtgctgaagctcgtatcaacaccagtggttggtgacctc  
c

ggcaccaccgaagccaacgcgtggaaatccaccctggttggtcacgacaccttcaccaaagttaaaccgtccgctgcttcccatcaccatcacca cca  
tt  
aataa

7/27/2015 – 7/31/2015

1. Retrieved the cellulose
2. Looked into ways to process cellulose into finer sheet
  - a. Tina used DIY paper-making protocol to do cellulose binding
3. Looked into ordering the sequence for two constructs:

CIPA: T7 promoter – RBS – CIPA CBD – Monomeric Streptavidin – End

BirA: J23104 promoter – RBS – AP Tag (Acceptor Peptide) – SgfI – aeBlue (Protein of Interest) – PmeI – RBS – biotin ligase – End

8/3/2015 – 8/7/2015

1. IDT synthesis arrived:
2. PCR cipA and birA
3. Digested and ligated into psb1C3 plasmid
4. Transformation of cipA and birA

8/10/2015 – 8/14/2015

1. Too many red colonies for cipA and birA: understandable because use RFP backbone
2. Picked out white colonies for overnight culture
3. cipA white colonies were actual gene; however, birA failed
4. Re-PCR and transformed of birA gene

8/17/2015 – 8/21/2015

1. Re-ligate and transformed with using CRATER method
  - a. Successful at removing most of the red colonies from using RFP as the backbone plasmid
  - b. However, the aeblue was not present
2. Colonies PCR and sequencing show that birA still is not there
3. Looked at geneious file for the gene sequence and realized that there is a PstI site at the birA gene.
4. Planned to use E and S as digestion enzyme instead of E and P

8/24/2015 – 8/28/2015

1. Digested with E and S digestion for birA; ligated and transformed
2. Sequencing did not work.
3. Couldn't figure out the problem

8/31/2015 – 9/4/2015

1. Tried to do digestion and ligation again but could not work
2. Looked into Gibson Assembly and tried Gibson
3. Sequencing after Gibson assembly still failed

9/7/2015 – 9/11/2015

1. Retried Gibson Assembly but still failed
2. Also did autolysis cloning and the autolysis cloning works successfully!

9/14/2015 – 9/18/2015

1. PCR the IDT Bira Gibson and send for sequencing
2. Sequencing of the PCR failed
3. This explained the reason that sequencing keeps failing
4. Did preliminary protein extraction of birA and cipA.
  - a. There are bands for cipA and birA; however, it does not match the right size of cipA and birA
  - b. Other protein extraction done at the same time failed. Possible reason for failure of cipA is that the flag-tag affinity does not work
5. Will redo protein extraction of cipA and re-clone birA