

# WESTER BLOT & BRADFORD ASSAY PROTOCOL

## SOLUTIONS

### Lysis Buffer (for whole bacteria lysis)

65mM Tris pH6.8 (7.4) → pH depends on the stacking gel

2% SDS

10% Glycerol

### Loading Buffer

65mM Tris-HCl pH6.8 (7.4) → pH depends on the stacking gel

2% SDS

10% Glycerol

5% mercaptoethanol (add under the HOOD)

### 1X Tris-Glycine Running Buffer

25mM Tris

190mM glycine

0.1% SDS

Adjust to pH 8.3 (if necessary or 7.4 depending on the gel)

### 1X Transfer Buffer

25mM Tris

190mM glycine

20% Methanol

Adjust to pH 8.3 (if necessary); if protein is >80kDa then add 0.1% of SDS

### Porceau Staining Buffer

0.2% (w/v) Porceau S

5% glacial acetic acid

### 1X TBST Washing Buffer

20mM Tris

150mM NaCl

0.1% Tween-20

### Blocking Buffer

Milk (5%) in TBST      OR

BSA (3%) in TBST (if you will use a phospho-Antibody detection)

### Antibody Solution

Antibody is diluted in Blocking Buffer (e.g. 1:5000)

### Stripping Buffer

20mL 10% SDS

12.5 mL 0.5M Tris-HCl

67.5mL MilliQ water

0.8mL  $\beta$ -mercaptoethanol

## PROCEDURE

### **1. Whole Cell Lysis**

- Grow bacteria to an OD<sub>600</sub> about 0.5. Take 5 mL of the bacterial solution and centrifuge at 16,000 rpm for 5 min
- Remove the supernatant carefully
- Resuspend the pellet in 50  $\mu$ L of cold Lysis Buffer
- Boil the sample for 10 min at 95°C and freeze them at -20°C for 10min. Repeat this 3 times.
- Take out a small fraction to evaluate the protein concentration by performing a Bradford Assay

*Note:* After measuring the protein concentration you can aliquot the sample and store it at -20°C or -80°C for long-term storage

### **2. Perform Bradford Assay**

- Dilute the Bradford Reagent (Bio-Rad) in 1:5 dH<sub>2</sub>O (and filter it with Whatman 540 filter paper)
- Add 10-20 $\mu$ L of your protein extract to 1mL of the diluted Bradford Reagent
- Measure the blue colour at 595nm
- Prepare also a standard curve (0.1mg/mL-1,0 mg/mL) of a known protein (e.g. BSA)

*Note:* Rule of thumb is that 1mg/mL should have an OD<sub>595</sub> of about 1.

### **3. Protein Separation by Gel electrophoresis**

- Add 10  $\mu$ L of the sample to 10  $\mu$ L of Loading Buffer. Add to this 1  $\mu$ L of bromophenol blue.
- Boil the samples at 95°C for 5 min.
- Quick-spin the sample down
- Set-up the SDS-PAGE chamber by putting the pre-cast gel in the construct and then fill the tank with *Running Buffer*. Remove the comb in the tank.
- Load equal amounts (20  $\mu$ L) into the wells of a mini pre-cast gel (Biorad)
- Run the gel at 200V (max. mA) and let the blue marker run down shortly before it reaches the bottom (in the last third of the gel)- Approx 35-45min. Stop the run and get the gel out of it cast. (CAREFUL!!! Gel is very fragile! Don't break it!)
- Cut one corner to remember in which well you have loaded which sample.

#### **4. Transferring proteins from the gel to the membrane**

- Let the gel stand in the Transfer buffer for 15 min.
- Prepare the transfer sandwich that consists of 2x sponge, then filter paper, then gel, then nitrocellulose membrane (also cut), then filter paper, then 2x sponges. All steps are done IN A TRANSFER BUFFER BATH.

*Note:* Make clear to press out all bubbles and close the transfer cassette

- Put the transfer cassette in the transfer tank. Verify that the blot is on the cathode and the membrane directed to the anode.
- Place also a block in the tank.
- Transfer the proteins at 100V for 30min (this step can immensely influence the bands that you will see later). Alternatively you transfer at 10mA during overnight in the cold room.

#### **5. Antibody incubations**

- Briefly rinse the blot in water and stain it with Ponceau solution to check the transfer quality.
- Rinse off the Ponceau S stain with three washes with TBST.
- Block in 3% BSA in TBST (required for phospho-Antibodies) or 5% milk in TBST at room temperature for 1 hour.
- Wash three times with TBST.
- Cut the membrane into two pieces: one which will be incubated with the antibody targeting the target protein; one against a control protein.
- Incubate overnight in the primary antibody solution against the target protein at 4°C.

*Note:* As control you can use GAPDH-targeting antibody

*Note:* The antibody should be diluted in the blocking buffer according to the manufacturer's recommended ratio. At room temperature: primary antibody needs about about 1-3 hours depending on the antibody quality and performance.

- Recollect the antibody!
- Rinse the blot 3-5 times with TBST for 5 min.
- Incubate the secondary antibody (e.g. HRP-conjugated) for 1 hour at room temperature
- Recollect the antibody!
- Rinse the blot 3-5 times for 5 min with TBST

#### **6. Imaging & Data Analysis**

- Apply the chemoluminescent substrate to the blot according to the manufacturers recommendation
- Capture the signal!

- Use the program imageJ to read the band intensity of the target protein (in comparison to the loading control GAPDH)

#### **7. Stripping and reprobing (if the protein sizes are quite different not necessary!!!)**

- Warm the Stripping Buffer to 50°C
- Add the Stripping Buffer to the membrane in a container designated for stripping. Incubate at 50°C for up to 45min with some agitation
- Rinse the blot under running water for 1 hour
- Transfer the membrane to a clean container, wash 5 times for 5 min with TBST
- Block in 3% BSA or 5% Milk in TBST at room temperature for 1 hour.
- Incubate with the primary antibody overnight. (Recollect the antibody!)
- Rinse the blot for 3-5 times for 5 min with TBST.
- Incubate with secondary antibody (e.g. HRP-conjugated) for 1 hour at room temperature.
- Rinse the blot 3-5 times for 5 min with TBST.

#### **8. Imaging & Data Analysis (if the protein sizes are quite different not necessary!!!)**

- Apply the chemiluminescent substrate to the blot following the manufacturer's suggestions.
- Capture the chemiluminescent signals
- Use imageJ to read the band intensity of the loading control proteins.
- Use the loading proteins levels to normalize the target protein levels.