

## **yeaRp Experiment Protocol**

### **Growth Medium: Luria Broth (LB)**

The test samples were first grown in LB overnight at 37°C. They were then washed with 0.85% NaCl solution. Washed samples were then resuspend into different concentrations of medium in a 96-well deep well plate and were further grown at 37°C until the bacteria reached mid-log phase. The fluorescence output was then measured using an EnVision multilabel reader.

Filter used on EnVision Multilabel Reader:

- Absorbance: Photometric 595nm,
- Excitation: 485nm FITC,
- Emission: 535nm FITC,
- Mirror module: FITC (403) on top.

10µl of antibiotics was added to each medium.

### **Characterization of *yeaRp* dynamic range in LB**

The concentrations used for the characterization of *yeaRp* was from 0 to 50 mM nitrate, with intervals of 10 mM.

<b>Expected final nitrate concentration (mM)</b>	<b>Actual final nitrate concentration (mM)</b>	<b>LB (ml)</b>	<b>1M KNO<sub>3</sub> (µl)</b>
0	0	10	0
10	9.89	10	100
20	19.58	10	200
30	29.10	10	300
40	38.42	10	400
50	47.57	10	500

The concentration of the characterization of *yeaRp* was from 0 to 10 mM of nitrate, with intervals of 2 mM.

<b>Expected final nitrate concentration (mM)</b>	<b>Actual final nitrate concentration (mM)</b>	<b>LB (ml)</b>	<b>1M KNO<sub>3</sub> (µl)</b>
0	0	10	0
2	1.99	10	20
4	3.98	10	40
6	5.96	10	60
8	7.93	10	80
10	9.89	10	100

## Growth Medium: M9

### Characterization of the promoter dynamic range in M9 minimal medium

The concentrations used for the characterization of *yeaRp* was from 0 to 2000  $\mu\text{M}$  nitrate, with 10 folds increase for each interval.

Expected final nitrate concentration ( $\mu\text{M}$ )	Actual final nitrate concentration ( $\mu\text{M}$ )	M9 (ml)	1M $\text{KNO}_3$ ( $\mu\text{l}$ )
0	0	10	0
20	19.89	10	0.2
200	199.76	10	2
2000	1994.02	10	20

The concentrations of the characterization of *yeaRp* was from 0 to 500  $\mu\text{M}$  of nitrate, with intervals of 100  $\mu\text{M}$ .

Expected final nitrate concentration ( $\mu\text{M}$ )	Actual final nitrate concentration ( $\mu\text{M}$ )	M9 (ml)	1M $\text{KNO}_3$ ( $\mu\text{l}$ )
0	0	10	0
100	99.89	10	1
200	199.76	10	2
300	299.61	10	3
400	399.44	10	4
500	499.25	10	5

### Reference:

Sambrook, J., & Russell, D. W. (2001). Molecular cloning. A laboratory manual. Third. *Cold Spring Harbor Laboratory Press, New York*.