

## **Phosphate sensor Experiment Protocol**

We have done a characterization on *phoAp* using Luria broth (LB) medium. Quantitative characterization on the promoter was done by measuring the fluorescence signal intensity using an EnVision multilabel reader.

Filter used on EnVision Multilabel Reader:

- Absorbance: Photometric 595nm,
- Excitation: 485nm FITC,
- Emission: 535nm FITC,
- Mirror module: FITC (403) on top.

*E. coli* strain DH10B was used, and the concentration of the characterization of *PphoA* promoter was from 0 to 300  $\mu\text{M}$  phosphate, with an intervals of 50 $\mu\text{M}$ .

60 $\mu\text{l}$  of antibiotics was added to each medium.

### **Preparing test medium with different concentration of phosphate**

We have prepared a solution of M9 minimal medium (J. Sambrook & D.W. Russell, 2001) and a solution of M9 minimal medium with Tris replacing phosphate. Test medium with different concentration of phosphate (0, 10, 30, 50, 100, 150, 200, 250, 300  $\mu\text{M}$ ) were made by mixing the 2 solution in the following ratio.

<b>Final Phosphate concentration (<math>\mu\text{M}</math>)</b>	<b>M9 minimal medium (<math>\mu\text{l}</math>)</b>	<b>M9 minimal medium without phosphate (replaced by Tris) (ml)</b>
0	0.00	60
10	8.58	60
30	25.73	60
50	42.85	60
100	85.71	60
150	130.43	60
200	171.42	60
250	216.26	60
300	260.87	60

## Characterization

Bacteria hosting pSB1C3-BBa\_*phoAp*-I13504 (or pSB1C3-BBa\_*phoBR*-I13504), positive control and negative control were first grown overnight in 5 ml Luria broth (LB) medium containing chloramphenicol at 37°C. The bacteria were then washed twice with 3 ml M9 minimal medium without phosphate (replaced by Tris), containing ampicillin. Then, the cells were resuspended in 5ml M9 minimal medium without phosphate (replaced by Tris) to obtain a final OD<sub>600</sub> of 4.50.  $\mu$ l of the prepared cell suspension were then added into 950  $\mu$ l of test medium with different concentrations of phosphate (containing chloramphenicol) in the 96-well deep well plate and further incubate at 37°C until the OD<sub>600</sub> of the cells reaches the mid-log phase. The fluorescence output were then measured using EnVision multilabel reader.

The results were obtained by biological triplicates and technical triplicates.

## Reference

Sambrook, J., & Russell, D. W. (2001). Molecular cloning. A laboratory manual. Third. *Cold Spring Harbor Laboratory Press, New York.*