

Transcribing and translating mCherry without cells

1. Set up PCR for negative control, positive control (luciferase), and mCherry.
 - a. 2 microliter DNA template
 - b. 5 microliter Amino Acid mixture (mix equal volumes of two amino acid mixture minus methionine and amino acid mixture minus cysteine; mix gently before use)
 - c. 20 microliters S30 premix without amino acid (mix before use)
 - d. 15 microliters S30 extract, circular (mix gently before use)
 - e. total volume: 50 microliters (add nuclease free to get final volume)
2. Vortexed gently and centrifuged for 5s.
3. Incubated 37C for 1 h.
4. Stopped rxn by placing in ice bath for 5 minutes.
5. Added 100 microliter of luciferase in A1, A2, A3, A4.
 - a. Negative Control B1
 - b. mCherry C1
6. Analyzed results using Tecan infinite m1000 plate reader at 37 C.
 - a. Shaking 5 s, Amplitude 4 mm
 - b. Luminescence:
 - i. Attenuation: Automatic
 - ii. Integration time: 10000 ms
 - iii. Settle time: 1000 ms
 - iv. Label name: Luciferase
 - c. Fluorescence Intensity:
 - i. Excitation = 587 nm
 - ii. Emission = 610 nm
 - iii. Gain: Calculated from well C1 (mCherry)
 - iv. Flashes: Mode 1 (400 Hz) 50
 - v. Z-Position: Manual 20000 micrometer
 - vi. Label Name: mCherry
7. For mCherry 55 microliters dH2O to wells, 45 microliters of mixture + assay reagent
 - a. For luciferase: 50 microliters dH2O to wells,

8. Added 50 microliters of assay reagent to PCR tubes.
9. Added 100 microliters from PCR tube to well.
10. Used plate reader.