

Nick's Lab Notebook

iGEM 2015

6.22.15

To do

-Reconstitute primers

-PCR to create inserts with cutsites

Goals

Two PCR products using tdTomato as Template

1. LexA DBD + RFP-Transcription Factor

2. pFig2c + LexA DBD + RFP-Transcription Factor

> pFig2c is an alpha factor responsive promoter

>

> [LexA DBD + RFP-Transcription Factor] will be spliced into a plasmid backbone (pGEM22) which contains pAga1(alpha factor responsive promoter)

###[] Reconstituting primers

#126 >> LexA DBD FW (Xho1)

#127 >> LexA DBD RV (not1)

#128 >> pFig2c FW (Apa1)

###[] PCR of inserts "Fig" & "Lex"

Fig

0.5ul tdTomato (template)

2.5ul FW #128

2.5ul RV #127

1.5ul DMSO

18ul H2O

25ul phusion mm

Lex

0.5ul tdTomato (template)

2.5ul FW #126

2.5ul RV #127

1.5ul DMSO

18ul H2O

25ul Phusion mm

*Products ran on 1% agarose gel

> "fig" insert will contain pFig2c promoter ~ 2.7kb

> "lex" piece will have no promoter ~ 2kb

###[] Gel results

"fig" piece has correct band

"lex" piece not seen

6.23.15

###[] Digestion pFig2c part + plasmid backbone
pGEM22(plasmid)

10ul pGEM22
3ul cutsmart
16ul H2O
0.5ul Apa1
----- RT for 1 hr
0.5ul Not1
----- 37C for 1 hr

> pGEM22 contains pAga1

"Fig" (insert) >>> add to PCR product tube

5ul Cutsmart
0.5ul Apa1
----- Rt for 1 hr
0.5ul Not1
----- 37C for 1 hr

###[] New PCR "lex" piece from scratcg

"Lex"

0.5ul tdTomato (template)
2.5ul FW #126
2.5ul RV #127
1.5ul DMSO
18ul H2O
25ul Phusion mm

6.24.15

To do

-New PCR LexA DBD + RFP-TransCription Factor
-w/ & w/o DMSO

###[] PCR LexA DBD + tdTomato(RFP-TF)

w/ DMSO

0.5ul tdTomato (template)
2.5ul FW #126
2.5ul RV #127
1.5ul DMSO
18ul H2O
25ul Phusion mm

w/o DMSO

0.5ul tdTomato (template)
2.5ul FW #126
2.5ul RV #127
19.5ul H2O

25ul Phusion mm

* Products ran on 1% agarose gel

###[] Gel results

Visible bands in both PCR products at correct size
(~2kb)

###[] Digest PCR products

"Lex" and "Fig" parts gel extracted
Plasmid Digest (pGEM22) was not cut

6.25.15

> Hy130E used as a backup in case pGEM22 isn't cut
> both are pNH604 backbone, Hy130E is empty

###[] Plasmid Digests: Hy130E; pGEM22 ###

Hy130E "A" ----> HyA
10ul Plasmid
3ul Cutsmart
16ul H2O
0.5ul Apa1
----- Rt for 1 hr
0.5ul Not1
----- 37C for 1 hr

Hy130E "X"----> HyX
10ul Plasmid
3ul Cutsmart
16ul H2O
0.5ul Xho1
0.5ul Not1
----- 37C for 1 hr

pGEM22 "A"
10ul Plasmid
3ul Cutsmart
16ul H2O
0.5ul Apa1
----- Rt for 1 hr
0.5ul Not1
----- 37C for 1 hr

pGEM22 "X"
10ul Plasmid
3ul Cutsmart
16ul H2O
0.5ul Xho1
0.5ul Not1

----- 37C for 1 hr

###[] Gel Extraction Plasmid Digests ###
Both Hy130E strains cut correctly
pGEM22 not cut

###[] Ligation: [HyX + LexA DBD-RFP-TF] & [HyA + LexA DBD-RFP-TF] ###
HyA + Fig
2ul HyA
8ul Fig
1ul Ligase
2ul T4 Buffer
7ul H2O

HyA + Fig (negative control)
2ul HyA
1ul Ligase
2ul T4 Ligase
15ul H2O

HyX + Lex
1.5ul HyX
3ul Lex
1ul Ligase
2ul T4 Buffer
12.5ul

HyX + Lex (Negative Control)
1.5ul HyX
1ul Ligase
2ul T4 Buffer
15.5ul

*Samples kept at 16C overnight

6.26.15

[] Transformation: HyX & HyA Ligations ###
HyA ----> 1x, 8x
HyX ----> 1x, 8x

Plated on LB/CARB, left at RT over the weekend

6.29.15

###[] Colony PCR: HyX & HyA ###
HyX MM
70ul GoTaq
7ul FW #126
7ul RV #127

21ul H2O

HyA MM

70ul GoTaq
7ul FW #128
7ul RV #127
21ul H2O

6.30.15

###[] Colony PCR ran on 1% agarose gel ###

Faint bands from HyX colonies

No bands from HyA colonies

HyX colonies: 1,2,4,5,6

---> made into liquid cultures

> liquid cultures were taped to the side of the shaker for some odd reason..

> lost colonies 1,4 liquid cultures

###[] Back up Ligation Re-do ###

HyA

2ul T4 buffer
2ul HyA
8ul Fig
1ul ligase
7ul H2O

HyA (negative control)

2ul HyA
2ul T4Buffer
1ul Ligase
15ul H2O

###[] Colony PCR: 8 new colonies picked from HyA plate ###

9x mm

90ul GoTaq
9ul FW #126
9ul RV #128
27ul H2O

*looking for 3kb PCR product

> sample ran on 1% agarose gel

> colony 8 contains band

7.1.15

Miniprep liquid cultures

Nanodrop results

- HyA colony 8: 139.6 ng/ul
- HyX colony 5: 474 ng/ul
- HyX Colony 6: 699.6 ng/ul
- HyX colony 2: 549.9 ng/ul

7/2/15

Colony PCR HyX 8x plate ###
9x mm

90ul GoTaq
9ul FW #126
9ul Rv #127
27ul H2O

Ran on 1% agarose gel

7/8/15

Gel Extraction new HyA & HyX Digests ###
Gel fell apart-->faulty agarose gel stock?

New Digests

Hy130E "A" ---> HyA
10ul Plasmid
3ul Cutsmart
16ul H2O
0.5ul Apa1

----- Rt for half day

0.5ul Not1

----- 37C for overnight

Hy130E "X"---> HyX
10ul Plasmid
3ul Cutsmart
16ul H2O
0.5ul Xho1
0.5ul Not1

----- 37C for overnight

7.9.15

Gel Extraction HyA & HyX Digestion

- HyX
 - cut correctly
- HyA
 - cut correctly

Ligation