



PROTOCOL TO PROMOTER ACTIVITY TEST UNDER CARBON STARVING

1. Make a pre-inoculation of DH5a :: puspA:GFP // pJ23101 of 15 mL containing 15 μ L of 35mg/mL chloramphenicol. (1 tube for each medium: LacBlue and LB). Incubate overnight at 37°C under 200 rpm.
2. 1mL of each medium should be collected to be used as test. Centrifuge at 4000rpm for 10 min. After, add 5 mL of medium LacBlue in the falcon previously containing LB and resuspend cells in ice bath carefully. Return cells to shaker for 4h at 37°C under 200rpm. Keep falcons with resuspended cells in new LB to comparison.
3. Recover aliquots of 100 μ L per hour. Centrifuge at 5000rpm for 5 min and freeze the pellet.
4. Proceed protein extraction, as suggested by our protocol, by physical lysis using ultrasound. As provided by protocol “**PROTOCOL TO PRODUCTION OF [6xHis] - GFP**” (Steps 5 to 8).
5. Make following SDS-PAGE:

Promoter J23101

M LacBlueON 0h 1h 2h 3h 4h Soluble Insoluble

Promoter UspA

M LacBlueON 0h 1h 2h 3h 4h Soluble Insoluble
6. Test fluorescence of samples containing J23101 (J) and UspA (U) at final of experiment in a Black Elisa plate in Viktor (Perkin Elmer) fluorimeter.

1	2	3	4	5	6	7	8	9	10	11	12
B1	B2	B3	J1	J2	J3	J4	U1	U2	U3	U4	X A (Medium LB)
B1	B2	B3	J1	J2	J3	J4	U1	U2	U3	U4	X B (Medium BLUE)

X	X	X	X	X	X	X	X	X	X	X	X	C
X	X	X	X	X	X	X	X	X	X	X	X	D
X	X	X	X	X	X	X	X	X	X	X	X	E
X	X	X	X	X	X	X	X	X	X	X	X	F
X	X	X	X	X	X	X	X	X	X	X	X	G
X	X	X	X	X	X	X	X	X	X	X	X	H

SOLUTIONS RECIPES AND BUFFERS

Lysis Buffer (1L):

Tris – 1,21g

NaH₂PO₄ – 6,89g

NaCl – 5,84g

Adjust pH to 8.0, bring volume to 1L. Sterilize by autoclaving at 121°C / 15 min.

Store at 4°C.

BLUE Medium (1L)

Bromothimol Blue – 0.08g

Ammonium Sulfate – 2g

Monobasic Potassium Phosphate – 13.6g

LB (Lysogenic Broth) – 1 mL of a solution 25g/1L

Cloramphenicol – 1 mL (35 mg/mL)

Sterilize by filtration at 0.22um seringe filter. Store at 4°C.