

Ligation Protocol

1. Add 2 μL of each restriction digest product to the reaction tube.
2. Add 2 μL of T4 DNA Ligase buffer to the tube.
3. Add 1 μL of ligase to the tube.
4. Add 11 μL of water to the tube.
5. Allow the reaction to proceed at room temperature for one hour.

Component	Volume
Part A Digest	2 μL
Part B Digest	2 μL
Linearized Plasmid Digest	2 μL
T4 DNA Ligase Buffer	2 μL
DNA Ligase	1 μL
dH ₂ O	11 μL
Total	20 μL

*****If this ligation protocol does not work, equimolar amounts of part A and part B were not added. Therefore, part A and part B volumes can be changed to obtain the correct ligation product. Total volume should be kept at 20 μL .**