

Protocol for InterLab Study

NJAU_China

Construction

Use BioBricks to construct devices.

Perform double digestion using SpeI & PstI for promoters, XbaI & PstI for GFP .

The system is as followed:

Buffer 4 ul

ddH₂O 12ul

Enzyme1 2ul

Enzyme2 2ul

DNA 20ul

Transform recombined plasmids into E. coli K-12 DH5-alpha.

Grow Cells for Measurement

Grow following devices and controls:

- Device1: J23101 + I13504 (B0034-E0040-B0015)
- Device2: J23106 + I13504 (B0034-E0040-B0015)
- Device3: J23117 + I13504 (B0034-E0040-B0015)
- Positive Control: BBa_I20270
- Negative Control: 1.E. coli K-12 DH5-alpha without any plasmid added 2.E. coli K-12 DH5-alpha with BBa_R0040 transferred

Followed the protocol provided by iGEM:

Streak out agar plates

LB Agar supplemented with 35ug/ml Chloramphenicol

Streak out 1 plate per device and control

Incubate plates until individual colonies are clearly visible at 37 C

Inoculate liquid culture

Use test tubes (15mm*150mm), oriente at an angle in the incubator, temperature at 37 C, shaking at 200 rpm for 16-18 hours. (Ronghua SHZ-82A Rotary Type Water Bath Thermostated Oscilllator)

3mL Luria Broth supplemented with 35ug/ml Chloramphenicol in each test tube.

Set up biological replicates in triplicate

Plate Reader Measurement

Obtain initial OD600 measurement

TECHCOPM UV2300

Set the spectrophotometer to read OD600, Take the measurement and record it

Dilute samples to an OD600 of 0.5

Calculate the dilution required for each sample, dilute each sample.

Re-measure samples on OD600

If OD600 is within 5% of 0.5, proceed

If OD600 is outside that range, recalculate dilution and remeasure until it's within 5%

Measure samples

TECAN infinite M200 with i-control 1.9

WHB Black Polystyrene Microplates

Set the reader to measure GFP: ex 485, em 528.

Add 100ul to a well for each sample.

Build standard curve by measuring sodium fluorescein with certain concentrations(500, 375, 250, 125, 50, 25, 10, 5, 0 ng/ml).

Measure.