

# iGEM 2015 Interlab Study

Organized by Jacob Beal, Traci Haddock, Markus Gershater, et al.

**Raytheon**  
BBN Technologies



SYNTHACE ■■



2 0 1 5  
G I A N T  
J A M B O R E E

# Why do an interlaboratory study?




## Measurement in Synthetic Biology:

- How precisely can the behavior of a part be characterized?
- How much do *de facto* protocols for measurement vary?
- What are the dominant causes of variation in measurement?

***Measurement is fundamental to everything we do.***


# What we asked teams to do

- Build three constitutive GFP constructs  
*Or arrange alternatives, if these weren't possible*
- Culture & measure fluorescence
- 3 biological replicates (Extra: x 3 technical rep.)





# Worldwide Participation!



## Participating Teams

Aalto-Helsinki	Birkbeck	CU Boulder	Glasgow	London Biohackspace	NJAU_China	Rock Ridge Virginia	TecCEM HS	Tufts	Waterloo
Aix-Marseille	BIT	Czech_Republic	Harvard_BioDesign	LZU	Northeastern	SCUT	Tec_Monterrey	UCL	William and Mary
Amoy	Boston University	Duke	Hong_Kong-CUHK	Marburg	NRP-UEA	SDU-Denmark	Tokyo Tech	UCLA	WLC-Milwaukee
ANU-Canberra	Brasil-USP	Edinburgh	HUST-China	METU_Turkey	NTNU-Trondheim	SPS-Singapore	Toronto	UC San Diego	WPI-Worcester
ATOMS-Turkiye	Cairo_Egypt	EPF_Lausanne	HZAU-China	Minnesota	NU_Kazakhstan	Stanford-Brown	Trento	UFMG_Brazil	
Austin_UTexas	Carnegie Mellon	ETH Zurich	IISER_Pune	MIT	OUC-China	Stockholm	TrinityCollegeDublin	UMaryland	
BHSF_Beijing	CityU_HK	Exeter University	KU Leuven	Nanjing_NFLS	Oxford	SYSU-Software	TU Delft	Utah_State	
Bielefeld	Cork	Freiburg	Leicester	Nankai	Paris-Saclay	SZMS_15_Shenzhen	TU Eindhoven	Vanderbilt	
BIOSINT Mexico	CSU_Fort_Collins	Gifu	Lethbridge	NEAU-China	Pasteur	TecCEM	Tuebingen	Vilnius-Lithuania	


# What did people use?

---


- Instruments:
  - Plate Reader: 56
  - Flow Cytometer: 23
  - Microscope: 7
  - Spectrofluorimeter: 5
  - Other: 2
- Strains:
  - DH5-alpha: 50
  - TOP10: 14
  - BL21: 9
  - Other: 20
    - XL1-Blue, JM109, 10G, RP437, MG1655, DH10-Beta, KRX

# Above and Beyond!

## Multiple strains and instruments!



## New Custom Equipment!




*We love how iGEMers surprise us!*

# Results: Measurement Precision


## J23101 (High) vs. J23106 (Med.)

iGEM 2015 Interlab: J23101 / J23106



## J23101 (High) vs. J23117 (Low)


iGEM 2015 Interlab: J23101 / J23117




***High measurements more precise than low***

# Results: 2015 vs. 2014

**2014: J23101 (High) vs. I20260 (Med)**



**2015: J23101 (High) vs. J23106 (Med)**




***More specific protocol → tighter distributions***




# Results: Cause of Variation?

## Strain




## Instrument



***Instrument matters, especially for low values.  
Strain, not so much.***

# Results: Intrateam Variation

iGEM 2015 Interlab: J23101 / J23106



***High instrument variation on single samples!***

# Time for Feedback!

*What worked? What didn't work?*



*How can it be bigger and better next year?*