




Transformation of GFP






Protein Production- Reminder







How is genetic engineering done?






Characteristic of Plasmids




GFP- Green Fluorescent Protein



Insertion of Fluorescent Genes



GFP Reporter Gene





GFP as a Reporter Gene





- A visible and defined signal
- Non-toxic protein
- Used for a wide range of research
- Used for research on gene location and expression regulation
- Used for research on protein location within a cell


Cancer and Blood Vessel Research



Nerve and Brain Cell Research



Transformation Regulation



Turning Genes on and Off




Gene expression types:

Constitutive: The gene is expressed at every stage

Induced: The gene is expressed as a reaction to an environmental signal

Silenced: No gene expression



The Lac Operon

Without Lactose



**Active
repressor**




Lac mRNA not transcribed

With Lactose



**Neutralized
repressor**




β -Galactosidase

Permease


Transacetylase

Transcription Regulation

Lactose Operon




Arabinose Operon




Gene Regulation

Natural Operon




Synthetic Operon





Bacterial DNA

Bacterial Cell




Genomic DNA




Plasmid DNA

Gene Expression in the pGLO plasmid




What is Transformation?

The insertion of a new gene into an organism



Bacterial Transformation



Transformation Types


Electroporation

An electric shock causes the membrane to be more permeable to DNA

Heat Shock

Calcium chloride and repeated heat and cold shocks enable DNA to enter the cell.

Transformation Steps




Transformation Steps

- Centrifugation of the cell culture and resuspension in a liquid containing calcium chloride
- Addition of the engineered pGLO plasmid
- Incubate on ice
- Heat shock at 42°C
- Return to ice
- Recovery in liquid LB
- Plate on selective medium

Why do we need these steps?

CaCl₂ Solution:

The Ca²⁺ ions bind negatively charged plasmid DNA and make it easier for it to permeate the membrane.




The Transformation Process

Heat:

Increases membrane permeability

Recovery in liquid LB:

Enables expression of the antibiotic resistance gene



LB Medium



Includes all the basic nutrient needed for bacteria:



Minerals

Salts


Amino Acids

Sugars

Plating Bacteria



Plating Bacteria



Which plate will fluoresce when containing our bacteria?

