


# RESEARCH LOGBOOK

## BACONTRACEPTION UNIVERSITAS INDONESIA

### Introduction

We divide our work into two works for wet laboratory. These are wet laboratory in one laboratory in Japan and two laboratory in Indonesia. The practitioner in Indonesia, especially work for constructing SboA and toggle switch in *Bacillus subtilis* wet lab are Jessica Luis, Faudina Fitra, Adi Nugroho Harlianto, and I Gusti Ngurah Raka. Otherwise, The Practitioner in Japan, especially work for constructing sboA expression in *Escherichia coli*, is Yudhistira Oktaviandi. Our purpose is expressing contraception in *Bacillus subtilis*, but we also do some work in *Escherichia coli* because we want to get more data to support our project (eg. the feasibility SboA in *Escherichia coli*). First, we will give short explanation how we divide works to make contraception. This explanation is given by mind map at below.



The different work in both laboratory (Japan and Indonesia) is absolutely affecting our logbook. That's why, we make two part of our logbook.

### Understanding Our Logbook

- We have been explained the protocol of contraception project which contain information about activity, general information ingredients and tools we use, so we


just give other information in this logbook, these are tittle of work (the reader can read our protocol and step of work from 'protocol page'), date of work, the practician who do this work (we usually use the keyword name of practician), the result of work, and more adjusment information (e.g the detailed amount of ingredients we use).


- The keyword name of practician :
  - Jessica Luis (J)
  - Adi Nugroho Harlianto (A)
  - I Gusti Ngurah Raka (R)
  - Faudina Nurilla Fitra (F)
  - Yudhistira Oktaviandie (Y)

**Part One : Insertion of SboA gene and SboA-AlbG gene into E coli DH5alpha using pET9a**


The practian is Yudhistira Oktaviandi (Y)

Date	No.	Activity
<b>24/06/2015</b> (Wednesday)	1.	Re-suspending Synthetic Dried DNA.
	2.	Preparing Agarose Agar 1%
<b>25/06/2015</b> (Thursday)	1.	Amplification of SboA & SboA-AlbG by PCR.
	2.	PCR Product Purification (SboA gene)
	3.	Checking PCR purified product by agarose electrophoresis (SboA gene)

			
26/06/2015 (Friday)	1.	Checking PCR purified product by aga//rose electrophoresis (SboA-AlbG gene)	


		
	2.	Agarose Gel DNA Purification for PCR product (SboA-AlbG gene)
<b>29/06/2015</b> <b>(Monday)</b>	1.	Culture E coli DH5alpha in LB broth
	2.	Checking purified PCR product concentration by NanoDrop (ng/ $\mu$ L) : Purified SboA gene Sample 1 : 20.8 Purified SboA gene Sample 2 : 15.9 Purified SboA-AlbG gene Sample 1 : 62.1 Purified SboA-AlbG gene Sample 2 : 115 <i>SboA gene has low concentration because using too much glass beads (5 <math>\mu</math>L for 25 <math>\mu</math>L PCR products)</i>
<b>30/06/2015</b> <b>(Tuesday)</b>	1.	Plasmid Isolation (pET9a) from cultured E coli DH5 $\alpha$
	2.	Checking isolated plasmid (pET9a) concentration by NanoDrop (ng/ $\mu$ L) : Purified pET9a Sample 1 : 272.1 Purified pET9a Sample 2 : 196.4

<b>01/07/2015</b> <b>(Wednesday)</b>	1.	Double Digest PCR-purified product and vector DNA : 100 ng CutSmart Buffer : 0.5 µL Restriction Enzyme : 2 x 0.2 µL (NdeI & BamHI) dH <sub>2</sub> O : 3.1 µL Incubate 2 hour at 37°C
	2.	Agarose gel DNA purification for digested DNA <i>No band appear when see using transilluminator/GelNinja because reaction volume very small</i>
<b>02/07/2015</b> <b>(Thursday)</b>	1.	Double Digest PCR-purified product and vector DNA : 900 ng CutSmart Buffer : 2 µL Restriction Enzyme : 2 x 0.5 µL (NdeI & BamHI) dH <sub>2</sub> O : until volume reaction 20 µL in total Incubate in 37°C static incubator for minimum overnight (16-20 hr)
<b>03/07/2015</b> <b>(Friday)</b>	1.	Agarose gel DNA purification for digested DNA <i>Band appears clearly in transilluminator/GelNinja</i>
	2.	Checking purified-digested DNA by agarose electrophoresis


		 <p>8000 6000 4000 200</p> <p>Vector      Insert (SboA)      Insert (SboA-AlbG)</p> <p><i>No band appear, there's something wrong when agarose gel DNA purification</i></p>
06/07/2015 (Monday)	1.	<p>Double Digest PCR-purified product and vector</p> <p>DNA : 500 ng</p> <p>CutSmart Buffer : 5 <math>\mu</math>L</p> <p>Restriction Enzyme : 2 x 0.25 <math>\mu</math>L (NdeI &amp; BamHI)</p> <p>dH<sub>2</sub>O : until volume reaction 20 <math>\mu</math>L in total</p> <p>Incubate in 37°C static incubator for 1 hour</p> <p><i>This composition was following NEB digestion protocol (Quantitative)</i></p>
	2.	<p>Agarose gel DNA purification for digested DNA</p> <p><i>Glass beads was decreased from 5 <math>\mu</math>L to 3 <math>\mu</math>L. Band appears clearly in transilluminator/GelNinja</i></p>
	3.	<p>Checking purified-digested DNA by agarose electrophoresis</p>



		<i>No band appear</i>
<b>07/07/2015</b> <b>(Tuesday)</b>	1.	Double Digest PCR-purified product and vector DNA : 500 ng  CutSmart Buffer : 5 µL  Restriction Enzyme : 2 x 0.25 µL (NdeI & BamHI)  dH <sub>2</sub> O : until volume reaction 20 µL in total Incubate in 37°C static incubator for 1 hour
	2.	Agarose gel DNA purification for digested DNA  <i>Band appears clearly in transilluminator/GelNinja</i>
	3.	Checking purified-digested DNA by agarose electrophoresis  <i>No band appear</i>
<b>08/07/2015</b> <b>(Wednesday)</b>	1.	Amplification of SboA & SboA-AlbG by PCR.
	2.	PCR Product Purification (SboA gene)
	3.	Agarose Gel DNA Purification for PCR product (SboA-AlbG gene)
	4.	Checking purified PCR product concentration by NanoDrop (ng/µL) : Purified SboA gene : 86.3 Purified SboA-AlbG gene : 109.5
<b>09/07/2015</b> <b>(Thursday)</b>	1.	Double Digest PCR-purified product and vector DNA : 1000 ng  CutSmart Buffer : 5 µL  Restriction Enzyme : 2 x 0.5 µL (NdeI & BamHI)  dH <sub>2</sub> O : until volume reaction 20 µL in total Incubate in 37°C static incubator for 1 hour  <i>This composition was following NEB digestion protocol (Quantitative) and the amount of digested DNA increased</i>
	2.	Agarose gel DNA purification for digested DNA  <i>Band appears clearly in transilluminator/GelNinja</i>
	3.	Checking purified-digested DNA by agarose electrophoresis  <i>No band appear even though amount of DNA had been increased</i>






		
10/07/2015 (Friday)	1.	<p>Double Digest PCR-purified product and vector</p> <p>DNA : 900 ng</p> <p>CutSmart Buffer : 2 <math>\mu</math>L</p> <p>Restriction Enzyme : 2 x 0.5 <math>\mu</math>L (NdeI &amp; BamHI)</p> <p>dH<sub>2</sub>O : until mixture volume 20 <math>\mu</math>L in total</p> <p>Incubate in 37°C static incubator for minimum overnight (16-20 hr)</p>
11/07/2015 (Saturday)	1.	Move digestion mixture into freezer
14/02/2015 (Tuesday)	1.	<p>Agarose gel DNA purification for digested DNA</p> <p><i>Band appears clearly in transilluminator/GelNinja</i></p>
	2.	Checking purified-digested DNA by agarose electrophoresis






		 <p><i>Drowning agar in EtBr for 30 minutes and increasing amount of DNA for checking by electrophoresis from 1 <math>\mu</math>L to 3 <math>\mu</math>L, finally band appears. So maybe it's not problem within digestion or purification but problem in checking by electrophoresis</i></p>
15/07/2015 (Wednesday)	1.	<p>Ligation using Takara Ligation Kit</p> <p>1. Insert (SboA gene) : 2 <math>\mu</math>L</p> <p>Vector : 0.5 <math>\mu</math>L</p> <p>Takara Ligation mix : 2.5 <math>\mu</math>L</p> <p>2. Insert (SboA-AlbG gene) : 2.2 <math>\mu</math>L</p> <p>Vector : 0.3 <math>\mu</math>L</p> <p>Takara Ligation mix : 2.5 <math>\mu</math>L</p> <p>Incubate at 16 °C for 30 minutes</p>
	2.	<p>Checking purified PCR product concentration by NanoDrop (ng/<math>\mu</math>L) :</p> <p>Purified digested SboA gene : 38.6</p> <p>Purified digested SboA-AlbG gene : 102</p> <p>Purified digested plasmid : 67</p>
	3.	Transformation
16/07/2015 (Thursday)	1.	Move plate from incubator and check for colonies.

		 <p><i>No colony grow; next ligation and transformation using T4 ligase (NEB products)</i></p>
21/07/2015 (Tuesday)	1.	Ligation using T4 ligase Insert: 3 $\mu$ L (SboA & SboA-AlbG in different tube) Vector: 1 $\mu$ L T4 ligase buffer: 2 $\mu$ L dH <sub>2</sub> O: 13 $\mu$ L T4 ligase: 1 $\mu$ L Incubate for 6 hours at 16 °C
	2.	Transformation
22/07/2015 (Wednesday)	1.	Move plate from incubator and check for colonies. 


		<i>No colony grow; composition of ligation mixture wasn't appropriate, next ligation will use quantitative method as suggested as NEB ligation protocol</i>
<b>23/07/2015</b> <b>(Thursday)</b>	1.	<p>Ligation using T4 ligase</p> <p>Insert: 6.5 ng (SboA gene)</p> <p>Vector: 67 ng</p> <p>T4 ligase buffer: 2 µL</p> <p>dH<sub>2</sub>O: until mixture volume 20 µL</p> <p>T4 ligase: 1 µL</p> <p>Insert: 351.75 ng (SboA – AlbG gene)</p> <p>Vector: 67 ng</p> <p>T4 ligase buffer: 2 µL</p> <p>dH<sub>2</sub>O: until mixture volume 20 µL</p> <p>T4 ligase: 1 µL</p> <p>Incubate for 16 - 20 hours at 16 °C</p>
	2.	Culture <i>E coli</i> DH5alpha in LB broth
<b>24/07/2015</b> <b>(Friday)</b>	1.	Transformation
	2.	Plasmid Isolation (pET9a) from cultured <i>E coli</i> DH5a
	3.	<p>Digest pET9a with 3 different treatment in 3 different plate :</p> <ol style="list-style-type: none"> <li>1. Digest using NdeI only with no insert</li> <li>2. Digest using BamHI only with no insert</li> <li>3. Digest using NdeI &amp; BamHI with no insert</li> </ol>
	4.	Agarose purification of digested pET9a
<b>25/07/2015</b> <b>(Saturday)</b>	1.	<p>Move plate from incubator and check for colonies.</p> 

		 <p><i>No colony grow</i></p>
<b>27/07/2015</b> (Monday)	1.	Ligation of digested pET9a
<b>28/07/2015</b> (Tuesday)	1.	Transformation
<b>29/07/2015</b> (Wednesday)	1.	<p>Move plate from incubator and check for colonies</p>  <p><i>2 plates have colonies inside which cut by NdeI only and BamHI only; for other plate which cut by NdeI &amp; BamHI, there's no colony. It could be conclude that nothing going wrong with digestion and ligation.</i></p>
<b>03/08/2015</b> (Monday)	1.	Amplification of SboA gene & SboA-AlbG gene
	2.	Purification of SboA gene & SboA-AlbG gene
	3.	<p>Checking purified PCR product concentration by NanoDrop (ng/μL) :</p> <p>Purified SboA gene: 119.8</p> <p>Purified SboA-AlbG gene: 100</p>
<b>04/08/2015</b> (Tuesday)	1.	<p>Double Digest PCR-purified product and vector</p> <p>DNA : 1000 ng</p>

		<p>CutSmart Buffer : 5 <math>\mu</math>L</p> <p>Restriction Enzyme : 2 x 0.5 <math>\mu</math>L (NdeI &amp; BamHI)</p> <p>dH<sub>2</sub>O : until volume reaction 20 <math>\mu</math>L in total</p> <p>Incubate in 37°C static incubator for 1 hour</p>
	2.	<p>Purification of digested SboA gene, SboA-AlbG gene &amp; pET9a</p> <p>Purified pET9a : 130.3 ng/ <math>\mu</math>L</p> <p><i>Only pET9a successfully purified</i></p>
	3.	<p>Ligation</p> <p>Insert (SboA gene) : 5.5 ng</p> <p>T4 Ligase Buffer : 2 <math>\mu</math>L</p> <p>T4 DNA Ligase : 1 <math>\mu</math>L</p> <p>Vector : 50 ng</p> <p>dH<sub>2</sub>O : 11.89 <math>\mu</math>L</p> <p>Insert (SboA-AlbG gene) : 262.5 ng</p> <p>T4 Ligase Buffer : 2 <math>\mu</math>L</p> <p>T4 DNA Ligase : 1 <math>\mu</math>L</p> <p>Vector : 50 ng</p> <p>dH<sub>2</sub>O : 11.89 <math>\mu</math>L</p> <p>Incubate for 16 - 20 hours at 16 °C</p>
<b>05/08/2015</b> <b>(Wednesday)</b>	1.	Move ligation product into -30 °C freezer
	2.	Transformation
<b>06/08/2015</b> <b>(Thursday)</b>	1.	Move plate from incubator and check for colonies
		

		 <p><i>No colony grow</i></p>
07/08/2015 (Friday)	1.	<p>Ligation</p> <p>Insert (SboA gene) : 5.5 ng</p> <p>T4 Ligase Buffer : 2 <math>\mu</math>L</p> <p>T4 DNA Ligase : 1 <math>\mu</math>L</p> <p>Vector : 50 ng</p> <p>dH<sub>2</sub>O : 11.89 <math>\mu</math>L</p> <p>Insert (SboA-AlbG gene) : 262.5 ng</p> <p>T4 Ligase Buffer : 2 <math>\mu</math>L</p> <p>T4 DNA Ligase : 1 <math>\mu</math>L</p> <p>Vector : 50 ng</p> <p>dH<sub>2</sub>O : 11.89 <math>\mu</math>L</p> <p>Incubate for 16 - 20 hours at 16 °C</p>
08/08/2015 (Saturday)	1.	Move ligation product into -30 °C freezer
10/08/2015 (Monday)	1.	<p>Transformation</p> <p><i>Decreasing heat shock temperature to 40 °C</i></p>
11/08/2015 (Tuesday)	1.	<p>Move plate from incubator and check for colonies</p> <p><i>9 colonies grow in SboA plate; 6 colonies grow in SboA-AlbG plate</i></p>
	2.	Replicate colonies
12/08/2015 (Wednesday)	1.	<p>Move replica to -20 °C freezer</p> 







Result: After transformation of SboA gene and SboA-AlbG gene using pET9a, there are some colonies grow, but it hasn't check yet.

Name : Yudisthira Oktaviandie

Project : Insertion of SboA gene and SboA-AlbG gene into shuttle vector PTRKH3-ldhGFP

<b>03/08/2015</b> <b>(Monday)</b>	1.	Amplification of SboA gene & SboA-AlbG gene
	2.	Purification of SboA gene & SboA-AlbG gene
	3.	Culture bacteria containing PTRKH3-ldhGFP
<b>04/08/2015</b> <b>(Tuesday)</b>	1.	Isolation of PTRKH3-ldhGFP
	2.	Double Digest PCR-purified product and vector DNA : 1000 ng  CutSmart Buffer : 5 µL  Restriction Enzyme : 2 x 0.5 µL (NdeI & BamHI)  dH <sub>2</sub> O : until volume reaction 20 µL in total  Incubate in 37°C static incubator for 1 hour
	3.	Purification of digested SboA gene, SboA-AlbG gene & PTRKH3-ldhGFP
	4.	Ligation  Insert (SboA gene) : 2.8 ng T4 Ligase Buffer : 2 µL T4 DNA Ligase : 1 µL Vector : 50 ng dH <sub>2</sub> O : until 20 µL  Insert (SboA-AlbG gene) : 131.25 ng T4 Ligase Buffer : 2 µL T4 DNA Ligase : 1 µL Vector : 50 ng dH <sub>2</sub> O : until 20 µL  Incubate for 16 - 20 hours at 16 °C
<b>05/08/2015</b> <b>(Wednesday)</b>	1.	Move ligation product into -30 °C freezer
	2.	Transformation
<b>06/08/2015</b> <b>(Thursday)</b>	1.	Move plate from incubator and check for colonies

		 <p><i>Bacteria grow as a smear not a dot (single colony)</i></p>
<b>07/08/2015 (Friday)</b>	1.	<p>Ligation</p> <p>Insert (SboA gene) : 2.8 ng</p> <p>T4 Ligase Buffer : 2 µL</p> <p>T4 DNA Ligase : 1 µL</p> <p>Vector : 50 ng</p> <p>dH<sub>2</sub>O : until 20 µL</p> <p>Insert (SboA-AlbG gene) : 131.25 ng</p> <p>T4 Ligase Buffer : 2 µL</p> <p>T4 DNA Ligase : 1 µL</p> <p>Vector : 50 ng</p> <p>dH<sub>2</sub>O : until 20 µL</p> <p>Incubate for 16 - 20 hours at 16 °C</p>
<b>08/08/2015 (Saturday)</b>	1.	Move ligation product into -30 °C freezer
<b>10/08/2015 (Monday)</b>	1.	Transformation
<b>11/08/2015</b>	1.	Move plate from incubator and check for colonies

(Tuesday)		 <p><i>Bacteria grow as a smear not a dot (single colony); there is some problem with MIC of antibiotics that we made. We're using erithromycin tablet instead of liquid erithromycin. Next time, we'll use liquid erithromycin to get adequate MIC of antibiotics.</i></p>
	2.	Isolation of plasmid from smear colonies
	3.	Check plasmid by digestion and run on agarose electrophoresis <i>From four plate we test, only plate 3 that contain plasmid.</i>

Result: Bacteria grow as a smear because of wrong concentration of antibiotics. From four plate, only one plate (plate 3) that contains plasmid and it hasn't check yet.

### **Part Two : Insertion of SboA gene and SboA-AlbG gene into E coli DH5alpha using pET9a**

Date	Number	Activity	Result	Detailed Information
06/05/2015 Wednesday	1.	Making Nutrient Broth and Nutrient Agar (J)		
	2.	Streaking Bacillus Subtillis		We incubate


		from Frozen Stock (F)		Bacillus subtilis from frozen stock at incubator 32°C We made 120 ml nutrient agar (2,76 grams) and 50 ml Nutrient Broth (0,4 grams)
<b>07/05/2015</b> <b>Thursday</b>	1.	Genom Isolation from Bacillus subtilis (F and A)		
	2.	Ethium bromide electrophoresis (F and A)		
<b>11/05/2015</b> <b>Monday</b>	1.	Making Mrs Broth and Mrs Agar (F)		
13/05/2015 Wednesday	1.	Isolation Lactobacillus Fermentum and Lactobacillus		


		Plantarum (R)		
	2.	Ethibium bromide electrophoresis Bromide for the result of isolation lactobacillus fermentum and plantarum.(R)		
<b>14/05/2015 Friday</b>	1.	Polimerase Chain Reaction with Taq Polymerase for L.plantarum, L.Fermentum, and Bacillus subtilis (A)		
	2.	Ethibium bromide electrophoresis for PCR confirmation (A)		
<b>19/05/2015 Wednesd ay</b>	1.	Doing PCR again for L.plantarum, L.Fermentum, and Bacillus subtilis. (F)		
	2.	Ethibium bromide electrophoresis for PCR confirmation(R)		
<b>6/07/2015 Monday</b>	1.	Inokulasikan bac A and Bac B in Nutrient Broth (R)		
	2.	Isolation plasmid		




		Bac A and Bac B (R)		
	3.	Ethidium bromide electrophoresis for Isolation confirmation (A)		
<b>07/07/2015 Tuesday</b>	1.	Streaking bacillus subtilis in Nutrient Agar and Lactobacillus fermentum and Lactobacillus plantarum in agar (J)		
	2.	Preparing materials for gateway (R)		
<b>08/07/2015 Wednesday</b>	1.	Gateway sboA (A)		
	2.	Primer gateway dilution (A)		
	3.	PCR gateway		
<b>09/07/2015 Thursday</b>	1.	Making Escherichia coli DH5alpha competent cell (J)		
	2.	Transformation BP into Escherichia coli DH5alpha (R)		
	3.	Incubate Escherichia coli containing BP		

		clonase (R)		
<b>10/07/2015 Friday</b>	1.	Making Escherichia coli DH5alpha competent cell (R)		
	2.	Transformation LR clonase into Escherichia coli DH5alpha competent cell (R)		
	3.	Incubate overnight the result of LR clonase transformation (R).		
<b>23/7/2015</b>	1.	Replica plating POXGW transformation in Escherichia coli DH5alpha (F)		

<p><b>24/07/2015</b></p>	<p>1.</p>	<p>Checking Replica plating POXGW (R)</p>		
<p><b>5/08/2015</b></p>	<p>1.</p>	<p>Transform PSB1C3 into Top 10 Escherichia coli (J)</p>		
	<p>2.</p>	<p>Streak bacillus subtilis overnight (A)</p>		

<b>11/08/2015</b> <b>Saturday</b>	1.	Isolation Escherichia coli containing BP clonase and electrophoresis isolation BP (J)		
	2.	PCR confiration using Pdonr (A)		

	3.	Isolation Escherichia coli containing LR clonase and electrophoresis LR clonase (R)		
	4.	Isolation Escherichia coli containing pOXGW (J)		
	5.	Replica plating for LR. BP, and POXGW. (J)		
<b>12/08/2015 Sunday</b>	1.	Isolation from all replica plating (J)		
<b>15/08/2015 Friday</b>	1.	Making Bacillus subtilis competent cell (R)		
	2.	Electroporasi bacillus subtilis		

PCR	sboA-taq
(R)	




Total  
Reaction  
: 20ul

**19/08/20**  
**15**  
**Wednesd**  
**ay**

Purification the  
result of PCR  
sboA-taq using  
Qiaex (J)

Digestion  
PSB1C3      linier  
and circular a (F)



**Digestio**  
**n**  
**Plasmid**  
**PSB1C3**  
**linier :**



4	ul
PSB1C3	
plasmid	
linier and	
4	ul
Enzym	
master	
mix	

Digestion  
POXGW (F)

Weight  
amount  
of  
POXGW  
: 184  
ng/ul, we  
use 5 ul



				(1 ug) BSA 2,5 ul, Buffer orange bioneer for digestion : 2,5 ul, XhoI and SpH I : 0,05 ul (each enzym restrictio n. Total reaction 25 ul
<b>20/08/2015 Thursday</b>	1.	Digestion SboA-Taq without heat kill (F)		
<b>22/08/2015 Saturday</b>	1.	Making E.coli competent cell (J))		
	2.	Heat Kill PSB1C3-SboA digestion on 20 /08/2015 (A)		
	3.	Transform PSB1C3-SboA into E.coli competent cell.		

		Then, incubate overnight (A)		
	4.	Isolation using High Speed Qiagen (J)		
<b>23/08/2015 Sunday</b>	1.	Check the PSB1C3-SboA transformation on 22/08/2015 (A)		
	2.	Replica plating PSB1C3-SboA transformation and Incubate overnight (A)		
<b>25/08/2015 Tuesday</b>	1.	PCR Accu Power Premix (R)		
<b>26/08/2015 Wednesday</b>	1.	Digestion PSB1C3-SboA using Pst 1 and EcoR1 (F)		
	2.	Nanodrop		


<b>31/08/2015</b> <b>Monday</b>	1.	Incubate PSB1C3-SboA (F)		
	2.	PCR confirmation PSB1C3-SboA (from isolation on 22/08/2015)		
	3.	Electrophoresis PCR confirmation PSB1C3-SboA sample (J)		
	4.	Digest PSB1C3-SboA (F)		
<b>01/09/2015</b> <b>Tuesday</b>	1.	Inoculate plasmid from replica plating into Luria Broth (F)		
	2.	Digestion using Cutsmart, but it hadn't be incubated in incubator 37°C (F)		Using 250 nanograms plasmid
<b>02/09/2015</b> <b>Wednesday</b>	1.	Electrophoresis the result of PCR and digestion on 31 August 2015 and on 01 September 2015 (A)		
<b>09/09/2015</b> <b>Wednesday</b>	1.	Drying DNA overnight		
	2.	Digestion POXGW (F)		

	3.	Purification plasmid (J)		
	4.	Electrophoresis PCR and digestion POXGW (A)		
<b>10/09/2015 Thursday</b>	1.	Plasmid containing POXGW-SboA sent to IGEM (A)		
	2.	Electrophoresis Purification plasmid on 09 September 2015 (A)		
<b>11/09/2015 Friday</b>	1.	Preparing material for induction (A)		
	2.	Preparing 12 plate filled Luria Broth and kanamycin for transformation (A)		
	3.	Making Luria Broth containing kanamycin		
	4.	TransformationPE T9A into BL21CP. Then, incubate overnight (A)		
	5.	Ligation toogle switch with front fragment		

		“POXGW” and fragment “PSB1C3”. Then, incubate at room temperature overnight. (A)		
<b>12/09/2015 Saturday</b>	1.	Checking the plate containing transformation BL21CP (R)		
	2.	Making BL21CP competent cell (R)		
	3.	Transformation BL21 again (A)		
	4.	Gibson Assembly between POXGW and PSB1C3 (A)		
<b>13/09/2015 Sunday</b>	1.	Incubate 3-4 hours BL21CP in incubation shaker for transformation.		
	2.	Transform gibbon into BL21CP (A)		
<b>14/09/2015 Monday</b>	1.	Purification SboA using Qiaquick (J)		
	2.	Double digestion Poxgw using XhoI and Sph 1 as enzym restriction (F)		
	3.	Double digestion pSB1C3 linear and circular using		

		EcoRI HF and PstI (F)		
	4.	Purification the result of digestion (POXGW and PSB1C3) (J)		
	5.	Agarose gel Electrophoresis for confirmation (J)		
<b>14/09/2015 Tuesday</b>	1.	Making Bacillus subtilis competent cell (R)		
	2.	Transformation into Bacillus subtilis competent cell (R)		
	3.	Electrophorasi Bacillus subtilis 168 (R)		
	4.	Agarose gel Electrophoresis again for confirmation the result of purification (yesterday)		
	5.	Re- Double digestion pSB1C3 linear and circular using EcoRI HF and PstI (F)		
	6.	Purification the result from double		



		digestion psB1C3 linear and circular (J)		
	7.	Agarose gel Electrophoresis for confirmation (J)		
<b>15/09/2015</b> <b>Wednesday</b>	1.	Ligation POXGW and SboA (F)		
	2.	Inoculate single coloni e.coli top 10 into luria broth then culture overnight/Prepari ng for competent cell on Saturday (F)		
<b>17/09/2015</b> <b>Friday</b>	1.	Transform the ligation POXGW and SboA into E.coli (Y)		
	2.	Replica plating		
<b>18/10/2015</b>	1.	Running SDS page (A)		

	2.	Isolation PSB1C3- Toogleswitch		
	3.	Well difussion assay		